

DONJOY®

**ELBOW
CCMI Mark III™**

DJO Global
Vista Manufacturing Center
3151 Scott Street
Vista, CA 92081-8553 USA
Tel: (800) 336-6569
www.djoglobal.com

Fax To: (800) 457-4221

Email To: defiance.help@djoglobal.com

BILL TO:

Account Name: _____
Account #: _____ P.O. #: _____
Contact: _____ Phone: _____ Ext. _____
Address: _____
City: _____ State: _____ Zip: _____

SHIPPING INFORMATION:

Same as "BILL TO"

Name: _____
Address: _____
City: _____ State: _____ Zip: _____
 Next Day Air 2nd Day Delivery
 3rd Day Delivery Ground Other _____

PRODUCT SELECTION: CCMI Mark III™

Frame: DEFIANCE® Football/Collision

Hinge:

Elbow (CI only)

Strapping Sequence:

CI

Measurement taken by:

Cast

CCMI Mark III

Recommended Football Options

- Chamois Liners/Pads
- Reinforced Carbon Composite Frame
- Metal Caps — Power Caps
- Square Strap Ends
- Player Last Name, First Initial


- School/Team Color: _____

FOR DJO USE ONLY:

Order #: _____
Brace #: _____

PATIENT INFORMATION:

Name: _____
Age: _____ Height: _____ Weight: _____
Measurement: Right Left Reverse
Measured By: _____
 New Brace Remeasurement/Repair Refurbish



COLORS:

- Black Gloss Seafoam Gloss True Red Gloss Bright Blue Gloss
- Black Matte Emerald Gloss Fuchsia Gloss Cotton Candy Gloss
- Smoke Matte Volt Gloss Crimson Gloss Olive Drab Green Matte
- White Gloss Bumble Gloss Grape Gloss Flat Dark Earth Matte
- Navy Gloss Safety Gloss Velvet Gloss Desert Tan Matte
- Powder Gloss Red, White & Blue Gloss

METALLIC GLOSS:

- Metallic Black Metallic White Gloss Metallic Forest Gloss
- Metallic Gunmetal Gloss Metallic Bright Blue Gloss Metallic Gold Gloss
- Metallic Silver Gloss Metallic Navy Gloss Metallic Tru Red Gloss

MULTIPLE COLORS:

Bicep Frame: _____
Forearm Frame: _____
Fade (matte colors excluded): Yes No

GRAPHICS (extra charge):

Stars & Stripes (non-metallic) Stars & Stripes (metallic)

FEATURES:

- Suede Liners/Pads
- Chamois Liners/Pads
- Swiveling Strap Tabs*
- Metal Caps - PowerCaps (recommended for contact sports)
- Reinforced Carbon Composite Frame
- Installed Extension Stop: 0° 10° 20° 30° 40° (circle one)
- Installed Flexion Stop: 45° 60° 75° 90° (circle one)
- Neoprene Suspension Strap Kit

ACCESSORIES (extra charge):

- High Strength Stainless Steel Gears Air Condyle Kit
- Full Extension Lockout Stop Floam Condyle Pads
- *Lycra® Suspension Sleeve Atrophy Pad Kit
- Neoprene Suspension Sleeve
- Sports Brace Cover
- Sili-Grip Strap Pad Kit
- Condyle Pad Kit


For further instructions, please refer to the CCMI Mark III™ instructional video and instructions card included with the instrument.

CAUTION!

FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY THE ORDER OF A LICENSED HEALTH CARE PROFESSIONAL

CUSTOM CONTOUR MEASURING INSTRUMENT - MARK III USE INSTRUCTIONS

NOTICE: For optimum long-term fit of the DEFIANCE® or DEFIANCE® III, measurements should be taken when there is no edema or atrophy present.



APPLICATION INFORMATION:

1. To unfold, detach and rotate swing arms of the CCMI Mark III™ device. **Figures A1 and A2.**
2. Have the patient stand with their back facing the practitioner and the shoulder of the affected arm slightly extended. The elbow should be relaxed and the palm should be facing forward. Do not allow the elbow to be flexed beyond 15 degrees. Apply device directly to bare arm. **Figure B.**
3. Turn the instrument over so that it reads opposite of the limb being measured. Pull apart the knee width caliper and position on patient's elbow. Arrows on condyle shells should align with the top of the olecranon process and midline with the humerus. **Figure C.**
4. For short thigh (biceps measurement) and calf (forearm measurement), extend tape measure marked in centimeters located proximal to elbow. **Figure D.**
5. To establish proper compression for elbow width measurement, support lateral side of instrument and press load plunger on medial condyle until triangles are aligned. **Figure E.**
6. Record elbow width on order form. **Figure F.**
7. Record lateral offset on order form. **Figure G.**
8. Record biceps and forearm circumference on order form. **Figure H.**



T 800.336.6569 D 760.727.1280 F 800.936.6569
2900 Lake Vista Drive | Lewisville, TX 75067 | U.S.A.
DJOGlobal.com/donjoy

Copyright © 2020 by DJO, LLC
MKT00-8822 Rev B

Individual results may vary. Neither DJO, LLC, nor any of its subsidiaries dispense medical advice. Please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you. Medical records, along with this form, may be required in order to support medical necessity. Copyright © 2020 by DJO, LLC